

Navigator 2.0 Accessories	
PRODUCT	PART NUMBER
Navigator 2.0 Roll Stand and Pole Clamp	N2-8800-01
Navigator 2.0 Storage & Carrying Case	N2-8000-07
Sterilization Tray, 2.6"(6.6cm)x12"(30.5cm) (wireless probe)	WP-9000-TR
Sterilization Tray, 6.5"(16.5)x18"(45.7cm) (fits all probe sizes and cable)	SP-9000-TR
Wireless Probe Battery Holder (1/BX)	WP-9050-01
Wireless Probe End Cap (1/BX)	WP-2000-10
Cable for CdTe Cabled Probes	GP-4001-00
Wireless Probe Cover (10/BX)	WP-9003-DRT
Cabled Probe Cover (10/BX)	SP-9004-DRC
Wireless Pilot Probe Batteries (12/BX)	WP-8500-12
Co-Pilot Sterile Device (12/BX)	GP-6801-00

Wireless Probe
Cover, Angled
Part # WP-9003-DRT



Sterilization Tray,
2.6"(6.6cm)x12"(30.5cm)
Part # WP-9000-TR



Navigator 2.0 Roll Stand
and Pole Clamp
Part # N2-8800-01



Wireless Probe Battery
(12/BX)
Part # WP-8500-12



Custom designed
sterilization trays to
hold probes in place

DILON NAVIGATOR 2.0

Sterilizable Surgical Gamma Probe



- Certified Sterilization Procedures
- High Sensitivity for Tc99 and I125
- Easy to Use
- Completely Wireless System
- 10+ Hours Lithium Ion Battery

Advanced Surgical Probe Technology

Optimized to facilitate precise localization in a broad range of anatomy.
Certified sterilization procedures for superior infection control.

Simple design means instantly ready;
no warm up or calibration

Head and Neck

Precise localization of parathyroid adenomas requires a probe designed with the unique requirements in the specific anatomy.

11mm Head & Neck Probe

Pulmonary Tumors

Addresses the challenges of localizing small, difficult to visualize lesions in MIS pulmonary procedures. Eliminates guesswork on tissue margins and target location. Robotic system compatible.

Daniel VATS Probe

Advanced Lymphatic Mapping

Both wireless and wired probes outperform competitive probes. The internal collimator and detector design gives excellent performance in multiple applications.

14 mm Angled Probe

14 mm Straight Probe

Laparoscopic Probes

Commonly used in radio-guided minimally invasive surgical procedures for thoracic and abdominal malignancy.

10 mm Short Laparoscopic Probe

10 mm Long Laparoscopic Probe

CONSOLE

Lightweight and compact console can easily be moved anywhere and can also be attached to a dedicated roll stand

DISPLAY

High visibility LED and antiglare display makes it easy to see from anywhere in the room

RANGE

Audible range to customize the sound profile for the amount of activity present

Certified sterilization options: Steris, STERRAD, TSO3

PROBES

Specialized probes remove the need for background filtering

COLLIMATOR

Built in collimator improves performance and never gets lost

WIRELESS

Zigbee wireless technology reduces potential Bluetooth interference

PROBE BATTERY

Convenient off the shelf batteries are easy to source, easy to change

High sensitivity in count detection



Works with BOTH Tc99 and I125

BATTERY

10+ hours of battery life; no cords to trip over

POWER BUTTON

Instant Power On Button



PROBE BATTERY

Convenient off the shelf batteries are easy to source, easy to change